



**Title:** *Hybridoma Generation: Cell Fusion of myeloma SP2/0*

No: RTLP-GL-Ab-18

Location:  
*Old CCRC Tripp Lab*

Approval Date:  
10 September 2004

Supersedes Date:

**Materials:**

Page 1 of 3

- |                        |                               |                           |                  |
|------------------------|-------------------------------|---------------------------|------------------|
| •Lab coat              | •37°C H <sub>2</sub> O bath   | •50 mL tubes              | •Pipettes        |
| •Gloves                | •Spleens from boosted animals | •Tissue Culture flasks    | •Pipetteman      |
| •RTLTP # for Media     | •Centrifuge                   | •37°C Incubator           | •Pipette Aid     |
| •RTLTP # for Fusion    | •Light microscope             | •96-well microtiter plate | •Pipetteman tips |
| •RTLTP # for ELISA     | •Hemocytometer                | •Syringe plunger          |                  |
| •Prepared SP2/0 cells* | •Trypan blue                  | •Mesh filter              |                  |

**\*1 week prior!**

**Procedure:**

**I. Cell Fusion and Hybridoma Selection:**

1. **One week prior to fusion** prepare myeloma SP2/0 cells in DMEM-10%/HEPES/pyruvate, e.g. day 14 prior to final boost. At day of fusion, you will need  $\sim 10^8$  SP2/0 cells/spleen to be fused
2. One day prior to fusion prepare fusion reagents:
  - a) DMEM-10% containing 10mM HEPES +1 mM sodium pyruvate (cDMEM-10)
  - b) DMEM-20% containing 10mM HEPES + 1 mM sodium pyruvate (cDMEM-20)
  - c) Serum-free DMEM (SF-DMEM)
  - d) SF-DMEM containing 10mM HEPES, and 1 mM sodium pyruvate (SF-cDMEM)
  - e) 50% PEG, sterile
  - f) ammonium chloride solution
  - g) cDMEM-20 + 1x HAT
  - h) cDMEM-20 + 1x HT

3. One day prior to fusion, split SP2/0 cells for vigorous growth

**4. Day of fusion (per spleen):**

- a) Check SP2/0 cells for contamination and growth.
- b) Pre-warm SF-DMEM, SF-cDMEM, cDMEM-20, and 50% PEG in a 37°C water bath.
- c) Harvest spleens in SF-DMEM.
- d) Dissociate spleens using mesh filter and wash in SF-DMEM.
- e) Lyse RBCs using ammonium chloride solution, e.g. resuspend spleen cell pellet in 5 ml of ammonium chloride solution, incubate 5 min at RT, wash 2 times in SF-DMEM
- f) Resuspend cell pellet in 10 ml of SF-cDMEM.
- g) Wash  $10^8$  SP2/0 in SF-DMEM 2 times.
- h) Resuspend cell pellet in 10 ml of SF-cDMEM.
- i) Using trypan blue, count viable spleen cells.
- j) Resuspend cells to  $2.5 \times 10^6$  cell/ml in cDMEM-20.

**6. Cell Fusion:**

- a) Mix SP2/0 cells and diluted spleen cells 1:1 in a 50 ml tube.
- b) Centrifuge to wash and pellet cells.
- c) Decant supernatant, and add 1 ml of pre-warmed PEG-50%, drop-by-drop over 1 minute, gently mixing the cells.
- d) Stir gently for an additional minute.
- e) Add 1 ml of pre-warmed SF-cDMEM to the cells, drop-by-drop over 1 minute, gently mixing the cells
- f) Repeat once with an additional 1 ml of pre-warmed SF-cDMEM.
- g) Add 7 ml of pre-warmed SF-cDMEM to tube, drop-by-drop with gentle stirring.
- h) At this point you should see cell clumping.
- i) Centrifuge lightly, 5 min at RT
- j) Decant the supernatant, and forcefully discharge 10 ml of pre-warmed cDMEM-20 into pellet to break up clumps

- k) Repeat “step j” several times to reach 40 mls of media, or roughly  $2.5 \times 10^6$  cell/ml calculated from “step 4.j” above.
- l) Place cells in a flask and incubate overnight at  $37^\circ\text{C}$  to allow fibroblast to adhere and be removed from the hybridomas.
- m) The next day, gently tap flask to dislodge hybridomas, and using a 10 ml pipet, add 2 drops/well of a 96-well plate using up all the volume of cells, and incubate overnight at  $37^\circ\text{C}$ .
- n) The next day, check seeded cells in wells. Wells should be nearly confluent.
- o) Add 2 drops of cDMEM-20 + HAT/well using a 10 ml pipet and incubate at  $37^\circ\text{C}$ .
- p) BE EXTREMELY CAREFUL to avoid contamination, using separate pipets for each 96-well plate, and don't mix up plate lids.
- q) On days 2-5, 7 and 9, aspirate half the volume of the plate well volume away, and replace with 2 drops/well cDMEM-20 + HAT from a 10 ml pipet – BE VERY ASEPTIC!
- r) Profound cell death should occur between day 2 and 3 when HAT is added. Between days 7 and 9, viable hybridomas should be detected growing. The feeding schedule is rigid for days 2-5 post-fusion!
- s) On day 14, repeat feeding except use cDMEM-20 + 1x HT.
- t) On day 15, repeat feeding using cDMEM-20 only, and screen the little buggers for antibody by ELISA by collecting supernatant.

Author	Management Approval/Date	Effective Date

