



**Title:** *Ammonium Sulfate Precipitation of Antibodies*

No: RTLP-GL-Ab-2

**Location:**  
*Old CCRC Tripp Lab*

**Approval Date:**  
10 September 2004

**Supersedes Date:**

**Materials:**

Page 1 of 1

- |                     |                              |                                    |   |
|---------------------|------------------------------|------------------------------------|---|
| • Lab coat          | • 4° C centrifuge            | • 50ml tubes                       | • Borate Buffer*  |
| • Gloves            | • 4° C cold storage          | • Saturated ammonium sulfate (SAS) |   |
| • Timer             | • distilled H <sub>2</sub> O | • Phosphate buffered saline (PBS)  | *Dependent upon whether product will be used as is or run on protein G column |
| • Aspirator         | • Microwave                  | • Dialysis tubing or cassettes     |   |
| • Transfer Pipette  | • Pipettes                   | • Spectrophotometer                |   |
| • 15ml tubes        | • Pipetteman                 |                                    |   |
| • 4° C centrifuge   | • Pipette aid                |                                    |   |
| • Ice bucket w/ ice | • Pipetteman tips            |                                    |   |
| • 4° C cold storage |                              |                                    |   |

**Procedure:**

1. Centrifuge monoclonal antibody supernatant for 30 minutes at 20,000 x g, 4°C. When centrifugation is complete, save the supernatant and store on ice.
2. Prepare saturated ammonium sulfate (SAS) by dissolving 760 grams of ammonium sulfate in one liter of boiling distilled water.
3. Let SAS cool and store at 4° C.
4. Keeping the SAS and supernatant on ice, add SAS drop-wise to the medium while stirring to 45% (vol/vol). Leave 2h at 4°C or overnight in cold room.
5. Centrifuge the mixture at 3000 x g for 30 minutes at 4°C, decant the supernatant and save precipitate (pellet).
6. Dissolve the precipitate in PBS or borate buffer\*.  
*NOTE: The buffer you choose depends on whether it will be used as is, or run through a protein G column.*
7. Dialyze against >20 volumes of buffer chosen for 24-48h at 4°C making sure to change the buffer at least 4 times during dialysis.
8. Determine the IgG concentration using A<sub>280</sub>.

