



**Title:** *Tissue Preparation Protocols*

No: RTLP-SARS-7

Location:  
*Old CCRC Tripp Lab*

Approval Date:  
16 August 2005

Supersedes Date:

**Materials:**

Page 1 of 2

•Lab coat	•Phosphate Buffered	•-80°C Storage	•Pipettes
•Gloves	Saline (PBS)	•dH <sub>2</sub> O	•Pipetteman
•Collection Tubes (K-EDTA)	•DMEM-5% FBS (if tissue prep)	•MicroChem Plus	•Pipette Aid
•Cryovials	•Pen/Strep (if tissue prep)	•Homogenizer (if tissue prep)	•Pipetteman tips
•Trizol	•L-Glutamine (if tissue prep)		•microfuge
			•microfuge tubes

### **Blood and Serum Collection for TCID<sub>50</sub> Assays and RNA Extraction**

1. Blood should be collected from animals in both K-EDTA (Purple top) and serum (Red top) collection tubes from the animals at the desired time points.
2. Whole blood is to be used for TCID<sub>50</sub> assays and RT-PCR. Whole blood (~0.5 ml) for TCID<sub>50</sub> assays should be placed into labeled cryovials and stored at -80°C. For RNA isolation, 100 µl of whole blood should be added to 900 µl of Trizol LS and the mixture stored at -80°C. ***Liberally spray the tubes with 10% Microchem Plus prior to removal from the BSL-3.***
3. Blood collected in serum tubes (Red top) should be allowed to clot for ~30 minutes, centrifuged at ~3,000 rpm for 10 minutes, and the serum removed and placed in a labeled cryovial for storage at -80°C.

### **Handling of Nasal Wash/Flush Samples for TCID<sub>50</sub> Assays and RNA Extraction**

1. Nasal washes/flushes will be performed on animals and the fluid collected at the time of necropsy.
2. Nasal wash fluid is to be used for TCID<sub>50</sub> assays and RT-PCR. The collected fluid (~0.15 ml) for TCID<sub>50</sub> assays should be placed into labeled cryovials and

stored at -80°C. For RNA isolation, 100 µl of nasal wash fluid should be added to 900 µl of Trizol LS and the mixture stored at -80°C. ***Liberally spray the tubes with 10% Microchem Plus prior to removal from the BSL-3.***

#### **Tissue Preparation (Lung, Kidney, Intestine) for TCID<sub>50</sub> Assays**

1. Frozen tissues should be sectioned, weighed, and 100 mg added to 1ml of DMEM containing phenol red, 5% FBS, P/S, L-glutamine. Tissue samples should then be homogenized with disposable grinders.
2. The homogenate should then be transferred to sterile microcentrifuge tubes, and centrifuged at ~2,000 rpm for 5 minutes.
3. Remove the homogenate fluid, place into 4 labeled cryovials, and store at -80°C.

#### **Tissue Preparation (Lung, Kidney, Intestine) for Total RNA Extraction**

1. Frozen tissues should be sectioned and weighed into 100 mg samples. Tissue samples should be homogenized using disposable grinders and ~1ml of Trizol reagent added to the tube (after homogenization). The sample volume should not exceed 10% of the volume of Trizol reagent used for homogenization.
2. Incubate the homogenized sample for 5 minutes at 15-30°C to permit complete dissociation of nucleotide complexes. Aliquot the sample into 4 labeled cryovials. Store at -80°C until transferred to Georgia Ali for RNA extraction. ***Liberally spray the tubes with 10% MicroChem Plus prior to removal from the BSL-3. Immediately transfer the samples to or store at -80°C.***